



Captair[®] Flex[™] S-M-SD-XLS



User Manual



Congratulations!

By choosing Captair[®] Flex[™] ductless filtering fume hoods, you have chosen an **efficient, reliable** way to ensure safety.

Your Captair[®] Flex[™] filtering fume hood guarantees that you are protected when working with chemicals that pose an inhalation risk for the user.

Your fume hood functions based on the recirculation of filtered air, which makes it possible to use high-efficiency filters to trap toxic particles and molecules and to recycle this air within the laboratory. The air downstream from the hood's modular filtration column, which uses the new Flex[™] technology, is thus free of all chemical pollutants.

Nevertheless, **the effectiveness of this system is directly dependent upon it being used correctly and monitored by its users.** Your laboratory may also benefit from ergonomic, economic, and ecological advantages provided by the Captair[®] Flex[™] fume hood throughout its life cycle.




**Do you want to work in a 100% safe environment year after year?
A careful reading of this manual is a must!**

Your new Captair® Flex™ ductless filtering fume hood offers several advantages from day one:


Protect the environment

Since no ductwork is required, a captair® Flex™ filtered enclosure totally eliminates the direct discharge of pollutants into the atmosphere and therefore contributes to the protection of the environment. Furthermore, a captair® Flex™ filtered enclosure does not generate any pollution linked to energy production unlike a traditional ducted fume hood.




Eliminate installation costs

The installation of a captair® Flex™ filtered enclosure is quick and easy. There is no need for a ductwork linked to an air supply / air extraction system, in comparison to traditional ducted systems. A single electrical outlet is all that is needed to make it work. Its setting up can be realized at any time, without any complex forward planning. Do not hesitate to compare this cost to a traditional ducted fume hood cost.




Achieve significant energy savings

Ductwork air balance is essential to the proper operation of a traditional ducted fume hood however; it is also the source of very important energy consumption. In fact, captair® Flex™ filtered enclosures do not generate any energy costs associated with the use of expensive extraction systems or conditioned air supplies. Operational costs remain at a minimum even when taking into account filter changes.



Benefit from the use of an immediately available unit, easy to relocate.

Captair® Flex™ filtered enclosures can be moved depending on the protection needs of the laboratory and can be easily relocated without disturbing room air balance.



Contents

Introduction	Page 2
---------------------	---------------

Getting started	Page 5
------------------------	---------------

- Flex[™] technology
- Description of the control box
- First start-up
- Navigating the digital display screen

Protecting yourself	Page 27
----------------------------	----------------

- Via the E.S.P. program (Erlab Safety Program)
- Via AFNOR NF X 15-211: 2009 standard

Maintenance	Page 35
--------------------	----------------

- Monitoring air face velocity
- Manually detecting filter saturation
- Automatically detecting filter saturation
- Replacing the filters
- The revolving system
- Cleaning and maintenance

Warning	Page 52
----------------	----------------

Getting started

- Flex™ technology
- Description of the control box
- First start-up
- Navigating the digital display screen

FLEX™ TECHNOLOGY

A modular filtration column which provides unprecedented adaptability and flexibility.

By combining molecular and particulate filtration technologies, erlab® was able to devise the new Flex™ modular filtration column technology.

This innovation developed by the erlab® R&D laboratory offers unprecedented flexibility, adaptability, and savings.












This single device can thus be configured to meet the protection needs of your laboratory.

Your device can be quickly reconfigured and can be easily used for other applications.

This flexibility was made possible by creating stackable, one-size-fits-all filtration cartridges—an innovation that is key to your new Captair® Flex™.


Possible configurations for your Captair® Flex™ fume hood

Manipulated chemicals/ Applications

	<p>Liquids</p> <p>For dilutions, titrations, extractions, transfers, ...</p>	<p>1C</p>  <p>Class 2 according to the NF X 15-211 : 2009 standard</p>	<p>2C</p>  <p>Class 1 according to the NF X 15-211 : 2009 standard</p>
	<p>Solids</p> <p>For sievings, grindings, weighings, formulations, compressions...</p>	<p>1P</p> 	
	<p>Liquids & solids</p> <p>For dissolutions, filtration, extractions, ...</p>	<p>1P 1C</p>  <p>Class 2 according to the NF X 15-211 : 2009 standard</p>	<p>1P 2C</p>  <p>Class 1 according to the NF X 15-211 : 2009 standard</p>
	<p>Handlings in clean rooms</p> <p>Up to class ISO 5 rooms</p>	<p>1C 1P</p>  <p>Class 2 according to the NF X 15-211 : 2009 standard</p>	<p>2C 1P</p>  <p>Class 1 according to the NF X 15-211 : 2009 standard</p>


You may contact your E.S.P. agent at any time to confirm the safety parameters related to the use of your device!

DESCRIPTION OF THE CONTROL BOX



① Master switch

② Model label




③ Name of the manufacturer or their representative



- 4 Valipass[®] form: Refer to the Valiquest[®] section for details.


- 5 Sampling port for the manual filter saturation test or for an ambient air sensor when the machine is equipped with the Molecode[™] S option.

- 6 Control panel




For information on the various scenarios related to the triggering of the device's indicator lights and audible alerts, refer to the **EVENT ALERTS** section of the **FIRST START-UP** chapter.

- 7 Fuse tag.



FIRST START-UP

You have carefully followed the steps described in the assembly instructions and your Captair® Flex™ fume hood is ready to use.




Turn on the device using the switch located on the left side of the control panel.
Activate the device's ventilation system using the ventilation button located underneath the digital display screen.

The green indicator light comes on and the screen displays the DEFAULT SCREEN DISPLAY.

The Timer™ begins to count down based on its factory settings (60 hours).

Default screen display



In order to guarantee that your device functions properly and in order to ensure your safety, your Captair® Flex™ fume hood is equipped with a monitoring system linked to the air flow parameters and to the device's Flex™ technology.

The Timer™

The Timer™ keeps track of the amount of time that your Captair® Flex™ fume hood has been in operation in order to remind the user when it is necessary to carry out a periodic filter saturation test.

The default setting (factory setting) is 60 hours (as required by NF X 15-211: 2009 standard).

An alarm will alert the user to check that the filter is still able to function properly.

A pre-warning alarm may go off between 0 and 10 hour(s) before the Timer™ alarm (according to factory settings).

Activation / Deactivation:

It is possible to activate or deactivate the Timer™ and its pre-warning alarm.

To make use of this functionality please refer to the section on “**Navigating the digital display screen.**”

Event alerts

This section describes all of the audible and visual alerts triggered by the alarms included in your Captair® Flex™ fume hood.

Event 1:

Fan control problem

Event conditions:

The fume hood does not reach the proper number of rotations per minute (RPMs) within a period of ten seconds after the ventilation system has been started.

Display screen



Alarm type



Yellow indicator light on, constant
Intermittent audible alarm

Solution:

Contact your maintenance service as quickly as possible.

Event 2:

Fan out of order

Event conditions:

The rotation speed of the fan is less than 700 RPM.

Display screen



Alarm type



Red indicator light on, constant
Continuous audible alarm

Solution :

Contact your maintenance service as quickly as possible.

Event 3:

Solvent detected in the sampling chamber

Event conditions:

The detection limit of the Molecode™ S unit has been reached.

Display screen



Alarm type



Yellow indicator light on, constant
Intermittent audible alarm

Solution :

Press the navigation and confirmation buttons simultaneously and release them to deactivate the alarm.

The red indicator light will stay on.

The concentration at the exhaust of a filter nearing its saturation point can build up very quickly. Consider replacing the filter.


Event 4:

Solvent detected 2

Event conditions:

The detection limit of the Molecode™ S unit has been surpassed.

Display screen



Alarm type



Red indicator light on, constant
Continuous audible alarm

Solution :

It is necessary to replace the filter.
Please contact your maintenance service.

Event 5:

Pollutant detected in the ambient air of the laboratory

Event conditions:

The ambient air sensor on the front detects that pollution is likely in the ambient air of the laboratory.

Display screen



Alarm type



Red indicator light on, blinking
Continuous audible alarm

Solution :

Press the navigation and confirmation buttons simultaneously and release them to deactivate the alarm.
Identify the source of the pollution.

Event 7:

Pre-warning alarm maintenance

Event conditions:

If the Timer™'s pre-warning alarm is activated, this alarm is triggered when the counter's value is less than the pre-warning value.

Display screen



Alarm type



Yellow indicator light on, blinking
Intermittent audible alarm

Solution :

Press the navigation and confirmation buttons simultaneously and release them to deactivate the alarm.

Consider checking whether or not the filter is saturated.

Refer to the maintenance section of this manual.


Event 8:

Maintenance

Event conditions:

The value of the Timer™ is equal to zero.

Display screen



Alarm type



Red indicator light on, constant
Continuous audible alarm

Solution :


Press the navigation and confirmation buttons simultaneously and release them to deactivate the alarm.

Consider checking whether or not the filter is saturated.


Refer to the maintenance section of this manual.

NAVIGATING THE DIGITAL DISPLAY SCREEN

DEFAULT SCREEN DISPLAY




ACCESSING THE MENUS



To access the menus of the screen, press down on the navigation and confirmation buttons simultaneously.

The menus will appear when the buttons are released.

I Main Menu



Notes:

Allows the user to view the status of the fume hood. Press the confirmation button to access the contents of the menu.

Allows the user to modify the parameters of the screen display language. Press the confirmation button to access the contents of the menu.

Allows the user to update information after testing for filter saturation. Press the confirmation button to access the contents of the menu.

Allows the user to modify the limit settings of the Molecode™ S detection sensor. Press the confirmation button to access the contents of the menu.


Allows the service engineer to set the ventilation flow rate. Protected by an access code. Press the confirmation button to access the contents of the menu.

Allows the user to activate or deactivate the Timer™. Press the confirmation button to access the contents of the menu.

Allows the service engineer to access the device's configuration. Protected by an access code. Press the confirmation button to access the contents of the menu.

Allows the user to exit the main menu. Press the confirmation button to display the default screen display.

I-1 Hood State



Notes :

This number corresponds to the number of fans detected by the device. It may be between 1 and 4 depending on the Captair® Flex™ model in question. (Refer to the technical specifications listed at the end of the manual.)

RPM: rotations / minute. This number indicates the rotation speed of each fan. It may be between 0 and 3000 depending on the Flex™ technology installed in your device and on the factory settings. The rotation speed must be the same for each fan installed, +/- 5%.

Indicates the presence or absence of the Molecode™ S option on your device.

This number indicates the sensitivity on the Molecode™ S detection sensor included on your device. If the device is not equipped with a Molecode™, the contents of this screen are not visible.

Display the Hood Stat menu.

I-2

Set Language

Notes:

By default the screen displays the name of the factory-set navigation language. Other possible display languages: French, Spanish, German, and Italian.

Display the current navigation language.

Display the Set Language menu.

Allows the user to modify the language displayed on the screen.


Confirm the chosen language to save the selection.

Confirm the chosen language to save the selection.


Confirm the chosen language to save the selection.

Confirm the chosen language to save the selection.

Display the Set Language menu.



I-3 Test filter



Notes:

N.B.: accessing this menu interrupts the ventilation system of your device. It is strongly recommended to access this menu after the filter saturation tests have been performed. Press the confirm button to proceed to the next step.

Parameter reset menu after a new filter has been installed.

Displays the month, year, and operating time of the device after a new filter has been installed in the Manufacturer menu.

Allows the user to confirm the installation date of the new filter and to reset the Timer™ if the Timer™ is activated according to factory settings. If the Timer™ is not activated, "H:Mn" is not displayed on the screen.

Allows the user to indicate the new filter installation date.

Displays the sensitivity value of the current sensor as stored in the device's memory or previously modified. If the Molecode™ S is not installed on the device, screens I-3-5, I-3-6, and I-3-7 do not display.

Allows the user to adjust the sensitivity value of the unit after a new filter has been installed. Scroll through the numbers using the scroll button and confirm each digit using the confirmation button.

Confirm the previously selected value.


Display the Test filter menu.

I-3 Test filter


I-3-2-1



I-3-2-2



I-3-2-3



Notes :


N.B.: Accessing this menu interrupts the ventilation system of your device. Press the confirm button to proceed to the next step.

After verifying the condition of the filter via a saturation test, display the installation date and operating time of this filter.

By pressing the confirm button, the Timer™ is reset according to the factory configuration if the Timer™ is activated.

Display the Test filter menu.

I-4 Set Threshold



Notes:


N.B.: accessing this menu interrupts the ventilation system of your device. Press the confirm button to proceed to the next step. **IMPORTANT:** this menu only appears if the device is equipped with a Molecode™ S.

Displays the sensitivity value of the Molecode™ S unit as configured in the factory or previously modified. Scroll through the numbers using the scroll button and confirm each digit using the confirmation button.

Allows the user to confirm the selected value.

Display the Set Threshold menu

I-5 Set Setpoint



Notes:

N.B.: the contents of this menu are protected by an access code. Contact your maintenance service for additional information.


By default, the screen displays 00000000. Using the scroll buttons, select the appropriate digit in the active box and press the confirm button to move to the next digit.

Confirm the selected access code.

Allows the user to enter the device's ventilation settings based on the Flex™ technology installed in your fume hood. This setting may be between 0 and 3000 RPM (rotations/minute). Using the scroll buttons, select the appropriate digit in the active box and press the confirm button to move to the next digit. Confirm the previously selected setting.

Display the Set Setpoint menu.

I-6 Set Timer



Notes:

This menu allows the user to activate or deactivate the Timer[™].

Displays the fact that the Timer[™] is activated.

Confirm activation of the Timer[™].


Display the **SetTimer** menu.

Displays the fact that the Timer[™] is deactivated.

Confirm deactivation of the Timer[™].

Display the **SetTimer** menu.

I-7 Manufacturer Menu



Notes:

N.B.: the contents of this menu are protected by an access code. Contact your maintenance service for additional information.

To enter the access code, scroll through the numbers in the active box using the scroll button. Press the confirm button to proceed to the next box. Once the code is entered, confirm the code using the OK button.

Indicates the number of fans detected by the CPU card when the device is first plugged in.


Indicates that the CPU card detects the Molecode™ S option when the device is first plugged in.

Counts the number of resets to the New filter mode (menu I-3-1-1).

Counts the number of resets to the Good filter mode (menu I-3-2-1).

Allows the user to access menu I-7-7.

I-7 Manufacturer Menu



Notes :

Allows the user to activate or deactivate the Timer™.


Allows the user to set the Timer™'s time value if the Timer™ is activated. Possible values: from 1 to 250 hours. Scroll through the numbers in the active box using the scroll button and confirm the selected value using the confirmation button.

Allows the user to activate or deactivate the pre-warning alarm for the filter saturation test.

Allows the user to set the point at which the pre-warning alarm is triggered. Possible values: between 0 and 10 hour(s). Scroll through the numbers in the active box using the scroll button and confirm the selected value using the confirmation button.

Allows the user to access menu I-7-11.

I-7 Manufacturer Menu



Notes:

Allows the user to modify the language displayed on screen.

Allows the user to modify the language displayed on screen.

Allows the user to modify the language displayed on screen.


Allows the user to modify the language displayed on screen.

Allows the user to modify the language displayed on screen.

Allows the user to access menu 1-7-12.

I-7 Manufacturer Menu

I-7-12



Notes:

Allows the service engineer to enter ventilation settings for the device. The rotation speed, between 0 and 3000 RPM (rotations / min), must correspond to the device's Flex[™] technology settings. Refer to the various Flex[™] technology settings for more information. Scroll through the numbers in the active box using the scroll button and confirm the selected value using the confirmation button.

The Debug ON mode allows the user to reduce the Timer[™] cycle in order to verify the proper functioning of the alarms in the factory (pre-warning, Timer[™], Molecode[™] S, etc.). The cycle takes 3 minutes to run.

Allows the service engineer to save or cancel parameters that were previously modified.

Display the default screen display mode.

Protecting yourself

- Via the E.S.P. program (Erlab Safety Program)
- Via AFNOR NF X 15-211: 2009 standard

THE **esp**® PROGRAM—KEEPING YOU SAFE

erlab® safety program

A long-lasting commitment by erlab® to the safety of the operator

An essential part of user safety, our laboratory's mission is to assess the interactive behavior of molecules and their interactive effect with the Flex™ filtration technology.

Based on this scientific investigation, your E.S.P.® specialist will recommend the most appropriate unit, define the adapted filtration column and the enclosure design and ensure a complete protection to the user. After installation, your E.S.P.® specialist will provide you with a constant monitoring of the unit regarding the handlings performed.

The **ValiQuest**® service:

Determine the most appropriate filtering fume hood for total safety during your handlings

Assisted by an E.S.P.® agent, you complete the investigation questionnaire, which precisely describes your intended chemical handlings. Our validation laboratory specialists will recommend the appropriate filtration fume hood and filter type. Personalized advice and accurate answer within 48 hours. A certificate validating the handling is supplied: real commitment of the manufacturer to the safety of the operator.

The **ValiPass**® service:


Certify and secure the usage framework at installation

When you receive your Captair® filtering fume hood, a usage certificate will give precise details on the chemicals to be used, the filter type and an estimation of its lifetime expectancy, for which your Captair® filtering fume hood has been validated. This certificate is a permanent reminder to the user or the safety officer of the data relating to their protection.

The **ValiGuard**® service:

A constant monitoring of your filtering fume hood.

Periodically (about every six months), the E.S.P.® agent will contact you to make sure that you have not changed your handlings and that the filter is still active. The E.S.P.® agent will show you how to perform step by step filter saturation tests and also the procedure for filter replacement. During this contact, if the E.S.P.® agent finds that there is a change in chemical handlings, you will be asked to complete a new questionnaire (see step 1). After review, a new certificate for use (to be placed on the front of the filtering fume hood) naming the approved chemicals will be sent to you to ensure that your chemical handlings are still performed within optimum safety conditions.



**Call your ESP® specialist today
and configure this unique Captair® Flex™
solution to your requirements.**

www.erlab.com

THE aliQuest® SERVICE

IDENTIFIES THE FILTERING FUME HOOD BEST SUITED FOR YOUR APPLICATIONS

Risk analysis and identification of protection needs

aliQuest® service

Tell us about your handlings.
Thank you to send us a completed file back to our sales department

A

Nr	Product name	Type of container	Open/Close	Duration (s)	Temp. (°C)	Handling frequency	Handling Quantity (ml ou g)	Duration (min)	Type of handling
1									
2									
3									
4	ABSOLUTE ALCOHOL	FLASK							
5	ACETALDEHYDE	FLASK							
6	ACETIC ACID	FLASK							
7	ACETIC ANHYDRIDE	FLASK							
8	ACETIC OXIDE	FLASK							
9	ACETONE	FLASK							
10	ACETONITRILE	FLASK							
11	ACETONE/ETHYLENE GLYCOL	FLASK							
12	ACETONE/ETHYLENE GLYCOL	FLASK							
13	ACETONE/ETHYLENE GLYCOL	FLASK							
14	ACETONE/ETHYLENE GLYCOL	FLASK							
15	ACETONE/ETHYLENE GLYCOL	FLASK							
16	ACETONE/ETHYLENE GLYCOL	FLASK							
17	ACETONE/ETHYLENE GLYCOL	FLASK							
18	ACETONE/ETHYLENE GLYCOL	FLASK							
19	ACETONE/ETHYLENE GLYCOL	FLASK							
20	ACETONE/ETHYLENE GLYCOL	FLASK							
21	ACETONE/ETHYLENE GLYCOL	FLASK							
22	ACETONE/ETHYLENE GLYCOL	FLASK							
23	ACETONE/ETHYLENE GLYCOL	FLASK							
24	ACETONE/ETHYLENE GLYCOL	FLASK							
25	ACETONE/ETHYLENE GLYCOL	FLASK							
26	ACETONE/ETHYLENE GLYCOL	FLASK							
27	ACETONE/ETHYLENE GLYCOL	FLASK							
28	ACETONE/ETHYLENE GLYCOL	FLASK							
29	ACETONE/ETHYLENE GLYCOL	FLASK							
30	ACETONE/ETHYLENE GLYCOL	FLASK							
31	ACETONE/ETHYLENE GLYCOL	FLASK							
32	ACETONE/ETHYLENE GLYCOL	FLASK							
33	ACETONE/ETHYLENE GLYCOL	FLASK							
34	ACETONE/ETHYLENE GLYCOL	FLASK							
35	ACETONE/ETHYLENE GLYCOL	FLASK							
36	ACETONE/ETHYLENE GLYCOL	FLASK							
37	ACETONE/ETHYLENE GLYCOL	FLASK							
38	ACETONE/ETHYLENE GLYCOL	FLASK							
39	ACETONE/ETHYLENE GLYCOL	FLASK							
40	ACETONE/ETHYLENE GLYCOL	FLASK							
41	ACETONE/ETHYLENE GLYCOL	FLASK							
42	ACETONE/ETHYLENE GLYCOL	FLASK							
43	ACETONE/ETHYLENE GLYCOL	FLASK							
44	ACETONE/ETHYLENE GLYCOL	FLASK							
45	ACETONE/ETHYLENE GLYCOL	FLASK							
46	ACETONE/ETHYLENE GLYCOL	FLASK							
47	ACETONE/ETHYLENE GLYCOL	FLASK							
48	ACETONE/ETHYLENE GLYCOL	FLASK							
49	ACETONE/ETHYLENE GLYCOL	FLASK							
50	ACETONE/ETHYLENE GLYCOL	FLASK							
51	ACETONE/ETHYLENE GLYCOL	FLASK							
52	ACETONE/ETHYLENE GLYCOL	FLASK							
53	ACETONE/ETHYLENE GLYCOL	FLASK							
54	ACETONE/ETHYLENE GLYCOL	FLASK							
55	ACETONE/ETHYLENE GLYCOL	FLASK							
56	ACETONE/ETHYLENE GLYCOL	FLASK							
57	ACETONE/ETHYLENE GLYCOL	FLASK							
58	ACETONE/ETHYLENE GLYCOL	FLASK							
59	ACETONE/ETHYLENE GLYCOL	FLASK							
60	ACETONE/ETHYLENE GLYCOL	FLASK							
61	ACETONE/ETHYLENE GLYCOL	FLASK							
62	ACETONE/ETHYLENE GLYCOL	FLASK							
63	ACETONE/ETHYLENE GLYCOL	FLASK							
64	ACETONE/ETHYLENE GLYCOL	FLASK							
65	ACETONE/ETHYLENE GLYCOL	FLASK							
66	ACETONE/ETHYLENE GLYCOL	FLASK							
67	ACETONE/ETHYLENE GLYCOL	FLASK							
68	ACETONE/ETHYLENE GLYCOL	FLASK							
69	ACETONE/ETHYLENE GLYCOL	FLASK							
70	ACETONE/ETHYLENE GLYCOL	FLASK							
71	ACETONE/ETHYLENE GLYCOL	FLASK							
72	ACETONE/ETHYLENE GLYCOL	FLASK							
73	ACETONE/ETHYLENE GLYCOL	FLASK							
74	ACETONE/ETHYLENE GLYCOL	FLASK							
75	ACETONE/ETHYLENE GLYCOL	FLASK							
76	ACETONE/ETHYLENE GLYCOL	FLASK							
77	ACETONE/ETHYLENE GLYCOL	FLASK							
78	ACETONE/ETHYLENE GLYCOL	FLASK							
79	ACETONE/ETHYLENE GLYCOL	FLASK							
80	ACETONE/ETHYLENE GLYCOL	FLASK							
81	ACETONE/ETHYLENE GLYCOL	FLASK							
82	ACETONE/ETHYLENE GLYCOL	FLASK							
83	ACETONE/ETHYLENE GLYCOL	FLASK							
84	ACETONE/ETHYLENE GLYCOL	FLASK							
85	ACETONE/ETHYLENE GLYCOL	FLASK							
86	ACETONE/ETHYLENE GLYCOL	FLASK							
87	ACETONE/ETHYLENE GLYCOL	FLASK							
88	ACETONE/ETHYLENE GLYCOL	FLASK							
89	ACETONE/ETHYLENE GLYCOL	FLASK							
90	ACETONE/ETHYLENE GLYCOL	FLASK							
91	ACETONE/ETHYLENE GLYCOL	FLASK							
92	ACETONE/ETHYLENE GLYCOL	FLASK							
93	ACETONE/ETHYLENE GLYCOL	FLASK							
94	ACETONE/ETHYLENE GLYCOL	FLASK							
95	ACETONE/ETHYLENE GLYCOL	FLASK							
96	ACETONE/ETHYLENE GLYCOL	FLASK							
97	ACETONE/ETHYLENE GLYCOL	FLASK							
98	ACETONE/ETHYLENE GLYCOL	FLASK							
99	ACETONE/ETHYLENE GLYCOL	FLASK							
100	ACETONE/ETHYLENE GLYCOL	FLASK							

B

Company : Post code :
 User's Name : Town :
 Address : Country :
 Tel : Fax :
 E-mail :

C

You require : ☒ a new ductless fume hood ☐ a replacement filter
 Please indicate the specific dimensions and/or the characteristics of the fume hood :
 Additional comments :

D

You require : ☐ a new ductless fume hood ☒ a replacement filter
 Model required :
 Serial number or inventory number :
 Additional comments :

ERLAB UK & Ireland
 34 Thomas's House
 34 Thomas's Square
 SP1 1BA Salisbury, Wiltshire
 Tel : 01722 341 940
 Fax : 01722 341 950
 E-mail : sales@erlab.co.uk

This service, which is completely free of charge, is offered to you before the purchase of your fume hood.

You may contact your E.S.P. agent at any time to reconfigure the Flex™ technology of your filtering fume hood. The fume hood can be modified according to the protection needs and the environment of your laboratory.

**This service is available online at www.erlab.com
or at one of our commercial branches.**

A Information related to your applications

- 1 Name of the chemical (drop-down menu): Select the name of the chemical handled by your lab from the list of options. If the name of the chemical does not appear on the list, select the «unlisted chemical» field from the list, then enter the name of the chemical. If this is the case, you must provide the material safety data sheet for the chemical in question to the erlab® laboratory so that we can evaluate the behavior of this molecule.
- 2 Type of container (drop-down menu): Select the type of recipient used to handle the chemical in question from the list of options.
- 3 Open / Closed (drop-down menu): Indicate whether the recipient used is open or closed.
- 4 Dilution (text field): Enter the concentration of the chemical in question as a weight percentage. If the chemical is used in its pure form, enter 100; if it is diluted, enter its concentration: 50, 1N, etc. and enter the name of the diluting agent on the next line.
- 5 Temperature (text field): Temperature of the chemical when handled (in °C).
- 6 Handling frequency (drop-down menu): Number of times that the chemical is actually handled per month. Select the frequency of use of the chemical in question from the list of options.
- 7 Handling quantity (drop-down menu): Volume of the chemical per use. Select the volume of the chemical in question from the list of options.
- 8 Duration (drop-down menu): Time, in minutes, per use. Select the length of time that the chemical in question is handled from the list of options.
- 9 Type of handlings (drop-down menu): Select the type of application performed from the list of options. If your application does not appear in the list, describe the application in the «Comments» text field of section C.

B Contact details : All fields marked with a * are required—this information allows us to properly process your file and to trace any information provided.

C Specifications of the request: New device, updating of your Flex™ filtration technology. The «Comments» text field allows you to provide all necessary information to enable us to refine the analysis conducted by the erlab® laboratory.


D Updating of your Flex™ filtration technology : If you select this option, you must indicate the brand name and type of fume hood that you currently own as well as its serial number.

ValiPass® certificate

The ValiPass certificate shall be delivered together with your Captairflex at the time of the first installation. It is the result of the ValiQuest you submitted to your supplier and is designed to remind the user(s) of this Captairflex unit of the chemical products which have been approved to be used, the expected life of the filter(s) and the saturation detection methods.

If this document is not delivered with your captairflex, we strongly recommend you to require it from your supplier, after filling the valiQuest questionnaire mentioned page 29.

If the application mentioned in the previous ValiQuest is designed to change, it is very important to fill in a new ValiQuest and require a new validation and a new ValiPass from your supplier before operating the change of application.



The image shows a yellow ValiPass certificate form with numbered callouts 1 through 9. The form includes fields for Captairflex type and Serial N°, checkboxes for handling types (Liquid chemicals, Liquid Chemicals + powder, Powders only, Clean room up to ISO 5), Validation Certificate N°, and dated. It also has sections for Chemical products authorized, Filter types, Installation date, Estimated life, Recommended chemical to test for a filter's saturation test with a reactive tube, and Molecode 5 detector installed? (Yes/Code, No).

- 1 Type of Captair® flex model:
- 2 Serial N° of the Captair® flex:
- 3 Main use of this Captairflex (depends on the Captair® flex model and the filters installed)
- 4 Restricted list of chemicals which can be used in this Captair® flex fume hood
- 5 Filter(s) type and date of installation
- 6 Estimated life of the filter (in months)
- 7 Chemical to be tested for the filter's saturation test with reactive tubes (check with suppliers such as Draeger, Gastec, Kitagawa, RAE system for the appropriate tube reference).
- 8 Inform if the Molecode automatic detector for the filter's is installed on this Captairflex and what is the setting Code N° (see page 40/41).
- 9 Alert the user to fill in a new ValiQuest questionnaire and transfer it to the supplier for a new validation, in case the application shall be changed, such as change of chemicals, quantities, etc.....

THE aliGuard® SERVICE :

CONTINUOUS MONITORING OF YOUR FILTERING FUME HOOD

The E.S.P.® agent contacts users about every six months.



The agent checks to see whether or not the laboratory's applications have changed.



aliQuest® Service

You are asked to describe your new application in order to validate the safety parameters.



The E.S.P.® agent will guide the user step by step through the tests to check for filter saturation as well as through the filter replacement procedure.

AFNOR NF X 15-211: 2009 standard

Appointed by the AFNOR, the Union de Normalisation de la Mécanique (UNM), composed of a college of experts (INRS, national organizations, professional syndicates), has established the AFNOR NF X 15-211 : 2009 standard. This standard applies to filtration fume hoods (also named recirculatory fume hoods or ETRAF) designed for research, analysis, education works,... for all laboratories where chemicals subjected to professional exposure limit values (OEL or TLV) are manipulated. This text imposes performances criteria linked to:

- **Filtration efficiency**
- **Containment efficiency**
- **Air face velocity**

As well as a specific documentation attached with each filtration fume hood.

The classes established by the standard


Class 1	Class 2
Filtration fume hood with safety reserve	Filtration fume hood without safety reserve
one main filtration level and one safety filtration level	One filtration level

Classification according to filtration type

	Denominations according to the NF X 15 211:2009 standard	Denominations equivalence for the erlab® products
Particles filtration*	Type P	Type P
Vapours filtration**	Type V	Type C
Particles and vapours filtration**	Type PV	Type PC

* :The particulate filter must be at least of H14 type according to the NF EN 1822-1 standard

** :The filters for vapours must be submitted to two successive performance tests with Cyclohexane and Isopropanol for the filters designed to retain Volatile Organic Compounds (VOC). Another test for the acid vapours is carried out with Hydrochloric acid.



The filtration efficiency

It is defined by the filter capacity to retain noxious molecules manipulated within the enclosure and qualifies the quality of the recirculated air at the filter exhaust.


	Class 1	Class 2
Normal operation phase	Emissions concentration at the filter exhaust must be lower than 1% of the TLV	
Detection phase	The concentration at the filter exhaust must be lower than 1% of the TLV, and the automatic saturation detector must warn the user	The concentration at the filter exhaust must be lower than 50% of the TLV
Safety operation phase	The concentration at the filter exhaust must be lower than 50% of the TLV ; its duration must not be lower than 1/12 of the normal operation phase duration.	

The retention capacities recorded during the tests performed on our filters demonstrate **the technical performance developed by ErLab®**.

These results ensure a very high protection level to the Captair® Flex™ users.

Example of test carried out on a Captair® Flex™ XL 714, equipped with BE+ filters in class 1

Isopropanol	Cyclohexane	HCL (35%)
2250 gr	3204 gr	7862 gr



Containment efficiency of the enclosure


It is defined by the fume hood capacity to maintain chemical vapors or particles within the enclosure without any propagation in the laboratory environment.

To prove this efficiency, a test is carried out following the protocol described in the standard.

SF₆ (Sulfur Hexafluoride) tracer gas emissions are made into the enclosure. A grid composed of sensors is placed in front of the working openings. Some samplings are done at the grid. Based on the emitted gas concentration and the samplings done (which will allow the determination

of an average operator exposure to this tracer gas), it is possible to determine a containment performance level of the filtration fume hood.

The containment threshold specified by the NFX 15-211:2009 standard imposes a maximum concentration of 0.1 ppm of SF₆ gas at the measuring points on the grid.



Air face velocity

It represents the fume hood capacity to create a dynamic barrier between the operator and the handling.

For filtration fume hoods with fixed front panel, the air face velocity at any point of the front openings must be between 0.4 and 0.6 m/s. They must be equipped with a permanent monitoring system for the air face velocity which also acts as a real time containment indicator.

The documentation

Filtration fume hoods must be delivered with a booklet containing an exhaustive list of chemicals, certified by the manufacturer, that can be handled into the filtering fume hood in the conditions described in the AFNOR NF X 15-211:2009 standard. The following information must be indicated in the booklet for each chemical listed:

- The chemical name, its formula, its CAS number, its boiling point, its molecular weight, its saturation vapour pressure
- The appropriate filter reference and its retention capacity during the normal operation phase
- The type of saturation detection system for the filter(s)
- The maximum quantity of the chemical that can be introduced within the enclosure
- The name of the test laboratory having performed the test



Erlab® has created its own booklet called the CHEMICAL LISTING. This booklet contains a list of approved chemicals indicating analysis data for about 700 molecular substances commonly used in laboratories. This booklet is delivered with each unit, as per the requirements of the AFNOR NF X 15-211:2009 standard.

Every Erlab® filtration fume hood is in compliance with these safety criteria.

The international standards

Erlab® products are in compliance with the following standards which guarantee your total protection.

France : **AFNOR NF X 15-211 : 2009**

USA : **ANSI/AIHA Z9.5
ASHRAE 110 : 1995**


Maintenance

- Monitoring air face velocity
- Manually detecting filter saturation
- Automatically detecting filter saturation
- Replacing the filters
- The revolving system
- Cleaning and maintenance

Monitoring air face velocity

The anemometer provided with your Captair[®] Flex[™] fume hood allows the user to constantly monitor the face velocity.

If the anemometer is not installed during assembly, please install this device according to the steps listed in the installation instructions provided with your fume hood.




When reading the face velocity, first be sure to check that the leveling bubble on the anemometer is positioned correctly. To ensure proper movement of the strip, the anemometer must be level. The air face velocity must be between 0.4 and 0.6 m/s (as required by AFNOR NF X 15-211: 2009 standard).

What should you do if the strip does not move or if its movement is irregular?

First of all, verify that the fume hood's ventilation system is turned on.

If the strip's movement is highly irregular, verify that the strip is correctly positioned in the housing located within the anemometer as shown in the diagrams below. The two pins of the strip must be correctly positioned in the support slots.



MANUALLY DETECTING FILTER SATURATION

Your application was analyzed and validated by the erlab[®] laboratory as part of the erlab Safety Program (E.S.P.). During this validation, the erlab[®] laboratory used the information provided to determine the chemical that must be the subject of the molecular filter saturation test. You have access to this information at any time on your Valipass[®] certificate.



Verify the chemical for which you will carry out the filter saturation test with a reactive tube



To carry out a manual filter saturation test, you need the following equipment:

- a pump (a simple plastic pump A is provided, For more accurate result use a professional pump type B from Draeger, Gastec, Kitagawa, RAE, etc...)
- a flexible silicone tube (provided)
- the chemical to be tested (not provided)
- an appropriate color coded tube for the chemical (Dreager, Gastec, Kitagawa, RAE, etc...)



Turn on the Captair[®] Flex[™] fume hood



Turn on the ventilation



Soak a roll up piece of paper in a container placed inside the enclosure. Pour a sufficient quantity of the chemical to be tested on the paper.



Connect the quick connector to the sampling port on the Control panel



Break off the 2 ends of the tube by inserting them into the plastic pump end.



Break off the 2 ends of the tube by inserting them into the designed hole of the professional pump



Connect the tube to the plastic pump following the direction of the arrow printed on the tube



Connect the tube to the professional pump following the direction of the arrow printed on the tube



Connect the other end of the tube to the silicone tube, pump gently to sample air from downstream. Refer to the instructions provided with the tube for the number of pump strokes



If you notice a change of coloration of the tube, your filter has to be changed.

After performing this test and depending on the results, you can then record this information by accessing the menus of the digital screen:

Either to enter the parameters related to the new filter installed (menu 1-3-1).



Or to confirm that the filter has not been changed (menu 1-3-2-1).



AUTOMATICALLY DETECTING FILTER SATURATION

Molecode™ S:

The automatic Molecode™ S detector is automatically activated when the ventilation system of your Captair® Flex™ fume hood is started.

It uses a sensor designed to detect various organic vapors and gases, which, based on the parameters set for the sensor, make it possible to detect saturation of the main filter.

How do you modify the settings of the sensor ?

Before turning on your Captair® Flex™ fume hood for the first time, verify that the detector setting corresponds to the code provided by the erlab® laboratory after analyzing and validating your application. This configuration was factory-set using the information provided.

This information is available at any time on your Valipass® form, which is attached to the control panel of your Captair® Flex™ fume hood.



Find the setting for the Molecode™ S, automatic detection of filter saturation by solvents.

To enter this setting, refer to section I-5 of the chapter on “Navigating the digital display screen.”

MODIFYING THE LIMIT SETTINGS FOR THE UNIT

For the detection of a single chemical:

Refer to the limits listed on the next page.

For the detection of multiple chemicals:

Contact your E.S.P.® agent as soon as possible in order to determine the limit value corresponding to your application.

MODIFYING THE LIMIT SETTINGS FOR THE MOLECODE™ S UNIT

You can modify the detection limit of the Molecode™ S at any time.

The table at left lists the setting codes to be entered for the detection of a single chemical.

To convey this information to the detector, refer to menu 1-4 of the section on “Navigating the digital display screen” in this manual.



CHEMICAL NAME	FORMULA	CAS Number	Molecode™ setting
1, 4-Dioxane	C4H8O2	123-91-1	3500
1,1,1-Trichloroethane	C2H3Cl3	71-55-6	2700
1,1-Dichloroethane	C2H4Cl2	75-34-3	2800
1,2-Dibromoethane	C2H4Br2	106-93-4	2700
1,2-Dichloroethane	C2H4Cl2	107-06-2	2700
1,2-Ethanediol	C2H6O2	107-21-1	3100
1-Aminobutane	C4H9NH2	109-73-9	2800
1-Butanol	C4H10O	71-36-3	4000
1-Chloro butane	C4H9Cl	109-69-3	3000
1-Propanol	C3H8O	71-23-8	4100
2, 2'-Dichlorodiethyl ether	C4H8OCl2	111-44-4	2800
2, 4-Dimethyl pentane	C7H16	108-08-7	3600
2,4-Dimethyl-3-pentanone	C7H14O	565-80-0	2800
2,6-Dimethyl-4-heptanone	C9H18O	108-83-8	3500
2-Amino butane	C4H9NH2	13952-84-6	2800
2-Butanol	C4H10O	78-92-2	4200
2-Butanone	C4H8O	78-93-3	4000
2-Butenal	C4H6O	4170-30-3	3600
2-Butoxyethanol	C6H14O2	111-76-2	2700
2-Chloroacetaldehyde	C2H3OCl	107-20-0	2700
2-Chloroethanal	C2H3OCl	107-20-0	2700
2-Chloroethanol	C2H5OCl	107-07-3	3000
2-Chloroethyl alcohol	C2H5OCl	107-07-3	3000
2-Ethoxyethanol	C4H10O2	110-80-5	3500
2-Furylmethanol	C5H6O2	98-00-0	3100
2-Heptanone	C7H14O	110-43-0	4300
2-Hexanone	C6H12O	591-78-6	2800
2-Hydroxymethylfuran	C5H6O2	98-00-0	3100
2-Methyl-1,3-butadiene	C5H8	78-79-5	2800
2-Methyl-1-propanol	C4H10O	78-83-1	2800
2-Methylbutane	C5H12	78-78-4	3500
2-Methylpropyl acetate	C6H12O2	110-19-0	3000
2-Methylpropyl ester of acetic acid	C6H12O2	110-19-0	3000
2-Pentanone	C5H10O	107-87-9	4000
2-Pentanone	C5H10O	107-87-9	4000
2-Phenyl propane	C9H12	98-82-8	3000
2-Propanol	C3H8O	67-63-0	4100
2-Propanone	C3H6O	67-64-1	4200
2-Propen-1-ol	C3H6O	107-18-6	2800
2-Propenol	C3H6O	107-18-6	2800
2-Propyl acetate	C5H10O2	108-21-4	3000
3-Chloro-1-propene	C3H5Cl	107-05-1	2800
3-Pentanone	C5H10O	96-22-0	4500
4-Methyl 2-pentanone	C6H12O	108-10-1	2800
4-Methyl-2-pentanone	C6H12O	108-10-1	2800

CHEMICAL NAME	FORMULA	CAS Number	Molecode™ setting
4-tert-Butyl toluene	C11H16	98-51-1	2800
5-Methyl-3-Heptanone	C2H6O2	541-85-5	3500
Absolute alcohol	C2H6O	64-17-5	4500
Acetaldehyde	C2H4O	75-07-0	2800
Acetic acid	C2H4O2	64-19-7	2800
Acetone	C3H6O	67-64-1	4200
Acetonitrile	C2H3N	75-05-8	3000
a-Chlorotoluene	C7H7Cl	100-44-7	2700
Alcohol	C2H6O	64-17-5	4500
Aldehyde ethylique	C2H4O	75-07-0	2800
Allyl alcohol	C3H6O	107-18-6	2800
Allyl chloride	C3H5Cl	107-05-1	2800
Allylglycidylether	C6H10O2	106-92-3	2800
Allylic alcohol	C3H6O	107-18-6	2800
Aminocyclohexane	C6H11NH2	108-91-8	2700
Aminoethane	C2H7N	75-04-7	2700
Ammonia	NH3	7664-41-7	2600
Ammonium hydroxide sol	NH4OH	7664-41-7	2600
Amyl alcohol n	C5H12O	71-41-0	3200
Azine	C5H5N	110-86-1	2700
Benzene chloride	C6H5Cl	108-90-7	2700
Benzyl chloride	C7H7Cl	100-44-7	2700
beta-Methyl acrolein	C4H6O	4170-30-3	3600
beta-Methylpropyl ethanoate	C6H12O2	110-19-0	3000
Bicyclopentadiene	C10H12	77-73-6	2800
Bromoethane	C2H5Br	74-96-4	2700
Bromoform	CHBr3	75-25-2	2800
Butyl acrylate	C7H12O2	141-32-2	2800
Butyl alcohol	C4H10O	71-36-3	4000
Butyl alcohol sec	C4H10O	78-92-2	4200
Butyl alcohol ter	C4H10O	75-65-0	3000
Butyl carbinol	C5H12O	71-41-0	3200
Butyl Cellosolve®	C6H14O2	111-76-2	2700
Butyl ether	C8H18O	142-96-1	3500
Butyl glycol	C6H14O2	111-76-2	2700
Butyl lactate	C7H14O3	138-22-7	3000
Butyl metacrylate	C18H14O2	97-88-1	2800
Butyl vinyl ether	C6H12O	111-34-2	3500
Butyl vinyl ether	C6H12O	111-34-2	3500
Butylene hydrate	C4H10O	78-92-2	4200
BVE	C6H12O	111-34-2	3500
Carbon disulfide	CS2	75-15-0	2600
Cellosolve®	C4H10O2	110-80-5	3500
Chlorobenzene	C6H5Cl	108-90-7	2700
Chlorothene	C2H3Cl3	71-55-6	2700

CHEMICAL NAME	FORMULA	CAS Number	Molecule™ setting
Cinamene	C8H8	100-42-5	2800
Crotonaldehyde	C4H6O	4170-30-3	3600
Cumene	C9H12	98-82-8	3000
Cumol	C9H12	98-82-8	3000
Cyanomethane	C2H3N	75-05-8	3000
Cyclohexane	C6H12	110-82-7	2800
Cyclohexanol	C6H12O	108-93-0	3000
Cyclohexanone	C6H10O	108-94-1	3500
Cyclohexene	C6H10	110-83-8	3800
Cyclohexyl alcohol	C6H12O	108-93-0	3000
Cyclohexyl ketone	C6H10O	108-94-1	3500
Cyclohexylamine	C6H11NH2	108-91-8	2700
Cyclopentane	C5H10	287-92-3	2800
Decane	C10H22	124-18-5	2700
Dibutyl ether	C8H18O	142-96-1	3500
Dichloropropane 1,2	C3H6Cl2	78-87-5	4500
Dicyclopentadiene	C10H12	77-73-6	2800
Diethamine	C4H11N	109-89-7	3200
Diethyl ether	C4H10O	60-29-7	3700
Diethyl ketone	C5H10O	96-22-0	4500
Diethyl oxide	C4H10O	60-29-7	3700
Diethylamine	C4H11N	109-89-7	3200
Diethylene dioxide	C4H8O2	123-91-1	3500
Diethylene oxide	C4H8O	109-99-9	4200
Diisopropyl ketone	C7H14O	565-80-0	2800
Diisopropyl ketone	C9H18O	108-83-8	3500
Dimethoxymethane	C3H8O2	109-87-5	3800
Dimethyl amine	C2H7N	124-40-3	2700
Dimethyl benzene (and isomers)	C8H10	95-47-6	2800
Dimethyl carbinol	C3H8O	67-63-0	4100
Dimethyl ketone	C3H6O	67-64-1	4200
Dimethyl sulfoxide	C2H6SO	67-68-5	2700
Dimethylacetone	C5H10O	96-22-0	4500
Dipropyl ketone	C7H14O	123-19-3	4300
DMA	C2H7N	124-40-3	2700
DMSO	C2H6SO	67-68-5	2700
Ethanal	C2H4O	75-07-0	2800
Ethanoic acid	C2H4O2	64-19-7	2800
Ethanol	C2H6O	64-17-5	4500
Ether	C4H10O	60-29-7	3700
Ethyl acetate	C4H8O2	141-78-6	4000
Ethyl acrylate	C5H8O2	140-88-5	3000
Ethyl alcohol	C2H6O	64-17-5	4500
Ethyl benzene	C8H10	100-41-4	2800
Ethyl bromide	C2H5Br	74-96-4	2700
Ethyl ethanoate	C4H8O2	141-78-6	4000
Ethyl ether	C4H10O	60-29-7	3700
Ethyl formate	C3H6O2	109-94-4	3000
Ethyl ketone	C5H10O	96-22-0	4500
Ethyl methyl ketone	C4H8O	78-93-3	4000
Ethyl nitrile	C2H3N	75-05-8	3000
Ethyl oxide	C4H10O	60-29-7	3700
Ethylamine	C2H7N	75-04-7	2700
Ethylen chlorhydrin	C2H5OCl	107-07-3	3000
Ethylene alcohol	C2H6O2	107-21-1	3100
Ethylene bromide	C2H4Br2	106-93-4	2700
Ethylene chloride	C2H4Cl2	107-06-2	2700
Ethylene chlorohydrin	C2H5OCl	107-07-3	3000
Ethylene dibromide	C2H4Br2	106-93-4	2700
Ethylene dichloride	C2H4Cl2	107-06-2	2700
Ethylene glycol	C2H6O2	107-21-1	3100
Ethylene glycol mono ethyl ether	C4H10O2	110-80-5	3500

CHEMICAL NAME	FORMULA	CAS Number	Molecule™ setting
Ethylidene chloride	C2H4Cl2	75-34-3	2800
Formic acid	CH2O2	64-18-6	2600
Furfuryl alcohol	C5H6O2	98-00-0	3100
Furyl carbinol	C5H6O2	98-00-0	3100
Glacial acetic acid (pure compound)	C2H4O2	64-19-7	2800
Glutaraldehyde	C5H8O2	111-30-8	2600
Glycol	C2H6O2	107-21-1	3100
Heptan-4-one	C7H14O	123-19-3	4300
Hexane	C6H14	110-54-3	2800
Hexone	C6H12O	108-10-1	2800
Hydroxybenzene	C6H6O	108-95-2	2800
Hydroxycyclohexane	C6H12O	108-93-0	3000
IPA	C3H8O	67-63-0	4100
Isoamyl acetate	C7H14O2	123-92-2	3700
Isoamyl alcohol	C5H12O	71-41-0	3200
Isobutanol	C4H10O	78-83-1	2800
Isobutenyl methyl ketone	C6H10O	141-79-7	3000
Isobutyl acetate	C6H12O2	110-19-0	3000
Isobutyl alcohol	C4H10O	78-83-1	2800
Isobutyl methyl carbinol	C6H14O	108-11-2	2700
Isobutyronone	C7H14O	565-80-0	2800
iso-Nitropropane	C3H7NO2	79-46-9	3000
Isooctane	C8H18	540-84-1	3400
Isopentane	C5H12	78-78-4	3500
Isopentyl acetate	C7H14O2	123-92-2	3700
Isoprene	C5H8	78-79-5	2800
Isopropanol	C3H8O	67-63-0	4100
Isopropyl acetate	C5H10O2	108-21-4	3000
Isopropyl alcohol	C3H8O	67-63-0	4100
Isopropyl benzene	C9H12	98-82-8	3000
Isopropyl benzene	C9H12	98-82-8	3000
Isopropylcarbinol	C4H10O	78-83-1	2800
Isopropylideneacetone	C6H10O	141-79-7	3000
MEK	C4H8O	78-93-3	4000
Mesityl oxide	C6H10O	141-79-7	3000
Mesitylene	C9H12	108-67-8	4000
Methanoic acid	CH2O2	64-18-6	2600
Methanol	CH4O	67-56-1	4200
Methyl acetate	C3H6O2	79-20-9	3000
Methyl acetone	C4H8O	78-93-3	4000
Methyl alcohol	CH4O	67-56-1	4200
Methyl benzene	C7H8	108-88-3	3000
Methyl butyl ketone	C6H12O	591-78-6	2800
Methyl cellosolve	C3H8O2	109-86-4	2700
Methyl chloroform	C2H3Cl3	71-55-6	2700
Methyl cyanide	C2H3N	75-05-8	3000
Methyl cyclohexanol	C7H14O	25639-42-3	3800
Methyl ethyl ketone	C4H8O	78-93-3	4000
Methyl formate	C2H4O2	107-31-3	2700
Methyl isobutenyl ketone	C6H10O	141-79-7	3000
Methyl isobutyl ketone	C6H12O	108-10-1	2800
Methyl metacrylate	C5H8O2	80-62-6	3000
Methyl propyl ketone	C5H10O	107-87-9	4000
Methyl styrene	C9H10	25013-15-4	3700
Methyl-2-propanol-2	C4H10O	75-65-0	3000
Methyl-3-butanol-1	C5H12O	71-41-0	3200
Methylal	C3H8O2	109-87-5	3800
methylamyl alcohol	C6H14O	108-11-2	2700
Methylethyl carbinol	C4H10O	78-92-2	4200
MIBC	C6H14O	108-11-2	2700
MIBK	C6H12O	108-10-1	2800
n-Amyl acetate	C7H14O2	123-92-2	3700


CHEMICAL NAME	FORMULA	CAS Number	Molecode™ setting
Naphta 30/60	85% Nonane/15% trimethylbenzene	8052-41-3	3600
n-Butanol	C ₄ H ₁₀ O	71-36-3	4000
n-Butyl acetate	C ₆ H ₁₂ O ₂	123-86-4	3000
n-Butyl amine	C ₄ H ₉ NH ₂	109-73-9	2800
n-Butyl chloride	C ₄ H ₉ Cl	109-69-3	3000
N-Ethylethanamine	C ₄ H ₁₁ N	109-89-7	3200
n-Hexane	C ₆ H ₁₄	110-54-3	2800
Nitroethane	C ₂ H ₅ NO ₂	79-24-3	3700
Nitropropane 2	C ₃ H ₇ NO ₂	79-46-9	3000
n-Octane	C ₈ H ₁₈	111-65-9	3400
Nonane all isomers	C ₉ H ₂₀	111-84-2	3600
n-Pentane	C ₅ H ₁₂	109-66-0	3000
o-Xylene	C ₈ H ₁₀	95-47-6	2800
p-Dichlorobenzene	C ₆ H ₄ Cl ₂	106-46-7	2700
Pentanol 1	C ₅ H ₁₂ O	71-41-0	3200
Phenol	C ₆ H ₆ O	108-95-2	2800
Phenyl chloride	C ₆ H ₅ Cl	108-90-7	2700
Phenyl hydroxide	C ₆ H ₆ O	108-95-2	2800
Phenylethane	C ₈ H ₁₀	100-41-4	2800
Phenylethylene	C ₈ H ₈	100-42-5	2800
Phenylmethane	C ₇ H ₈	108-88-3	3000
Propanol-1	C ₃ H ₈ O	71-23-8	4100
Propyl acetate	C ₅ H ₁₀ O ₂	109-60-4	3000
Propyl alcohol	C ₃ H ₈ O	71-23-8	4100
Propylacetate	C ₅ H ₁₀ O ₂	109-60-4	3000
Propylamine	C ₃ H ₉ N	107-10-8	3500
Propylbenzene	C ₉ H ₁₂	108-67-8	4000
Propylene aldehyde	C ₄ H ₆ O	4170-30-3	3600
Propylene chloride	C ₃ H ₆ Cl ₂	78-87-5	4500
Propylene dichloride	C ₃ H ₆ Cl ₂	78-87-5	4500
Propylene oxide	C ₃ H ₆ O	75-56-9	2800
Pyridine	C ₅ H ₅ N	110-86-1	2700
SBA	C ₄ H ₁₀ O	78-92-2	4200
sec-Amyl acetate	C ₇ H ₁₄ O ₂	123-92-2	3700
sec-Butyl amine	C ₄ H ₉ NH ₂	13952-84-6	2800
Stoddard solvent	85% Nonane/15% trimethylbenzene	8052-41-3	3600
Styrene	C ₈ H ₈	100-42-5	2800
ter-Butyl acetate	C ₆ H ₁₂ O ₂	540-88-5	3500
tert-Butyl alcohol	C ₄ H ₁₀ O	75-65-0	3000
Tetrahydrofuran	C ₄ H ₈ O	109-99-9	4200
THF	C ₄ H ₈ O	109-99-9	4200
TMA	C ₃ H ₉ N	75-50-3	2700
Toluene	C ₇ H ₈	108-88-3	3000
Toluol	C ₇ H ₈	108-88-3	3000
Tribromomethane	CHBr ₃	75-25-2	2800
Triethylamine	C ₆ H ₁₅ N	121-44-8	2700
Trimethyl carbinol	C ₄ H ₁₀ O	75-65-0	3000
Trimethyl pentane-2,2,4	C ₈ H ₁₈	540-84-1	3400
Trimethylamine	C ₃ H ₉ N	75-50-3	2700
Trimethylbenzene	C ₉ H ₁₂	108-67-8	4000
Turpentine oil	C ₁₀ H ₁₆	8006-64-2	3500
Vinyl acetate	C ₄ H ₆ O ₂	108-05-4	3500
Vinyl carbinol	C ₃ H ₆ O	107-18-6	2800
Vinyl toluene	C ₉ H ₁₀	25013-15-4	3700
Vinylbenzene	C ₈ H ₈	100-42-5	2800
White spirit	85% Nonane/15% trimethylbenzene	8052-41-3	3600
Xylene (isomers)	C ₈ H ₁₀	1330-20-7	2800

REPLACING THE FILTERS


Your device is equipped with FLEX™ technology that was configured to the user's protection needs when the device was purchased. The design of the column is thus directly dependent on the applications carried out in the enclosure. These applications may change over time. Your FLEX™ technology can therefore be reconfigured if your fume hood is used for applications other than those anticipated when the device is first set up.

The table below summarizes all possible Flex™ technology configurations for your device.


Molecular filter



HEPA filter H14



SD32I - SD 39I - SD 48I - SD 63I	IC	x1	
SD32I - SD 39I - SD 48I - SD 63I	2C	x2	
SD 32I	IP		x1
SD32I - SD 39I - SD 48I - SD 63I	IP IC	x1	x1
SD32I - SD 39I - SD 48I - SD 63I	IP 2C	x2	x1
SD32I - SD 39I - SD 48I - SD 63I	IC IP	x1	x1
SD32I - SD 39I - SD 48I - SD 63I	2C IP	x2	x1




Each molecular filter is labeled as shown at left. Please follow these markings.

The table below summarizes the different types of carbon filters that erlab® offers as well as the range of application of these filters.


Type AS	For organic vapors
Type BE +	For acid and organic vapors
Type K	For ammonia vapors
Type F	For formaldehyde vapors
HEPA H14	For powders




For these operations, we strongly recommend that the user or maintenance technician wear the necessary safety equipment, including: safety glasses, a safety suit, and gloves.



Turn off your Captair[®] Flex[™] fume hood.



Remove the two protective casings on either side of the hood.




You should be able to see the ventilation card located beside the master switch.




Disconnect the power supply to the ventilation card then the flexible silicone tube.

Remove all items comprising your fume hood's Flex[™] technology.


IC




IP




2C




IP IC




IC IP



IP 2C




2C IP




After carefully removing the filters from their packaging, assemble the column based on the Flex[™] technology to be installed.


IC




IP




2C




IP IC




IC IP




IP 2C



2C IP







The column is guaranteed to be airtight:

The design of the filtration cartridges guarantees that the filtration column is completely airtight. The vertical stacking principle of the modular filtration column uses a gel-based joint technology that, when applied to the perimeter of all of the filtration cartridges, ensures a perfect seal of the entire column thanks to gravity.


Be sure that the filters are positioned as shown in this diagram.




When reassembling the filtration column, be sure that the ventilation module is positioned correctly.



Reconnect the power supply to the ventilation card and the flexible silicone tube to the sampling chamber.



Re-attach the two protective casings.



You may then turn on your Captair® Flex[™] fume hood.

Do not forget to record the information related to the new filter installed by accessing the menus of the digital display screen (menu 1-3).




THE «REVOLVING» SYSTEM

This simple concept involves placing one chamber equipped with a molecular detection system between two filters having the same capacity. When the main filter is saturated, the molecules are directed to an identical back-up filter that is placed just above the detection chamber.

This system prevents all molecules from being released into the environment since they are automatically absorbed by the back-up filter. The back-up filter replaces the main filter when the main filter has reached its maximum saturation point. A new filter is then installed in place of the back-up filter. This cycle may be repeated indefinitely.

Compared to traditional filtration systems, the “revolving filter” concept increases the retention capacity of the filter by 25% and decreases replacement costs by 25%.



The “revolving system” may be used with following Flex[™] technology types:

IP 2C

2C

2C I P

CLEANING AND MAINTENANCE

CHECKING THE MECHANICAL PARTS

Hinges:

Hinges must be properly attached and should not jam; they must allow the front panel of the hood to be instantly and easily lifted upward.

Acrylic parts:

These parts must be clean; white streaks or spatters indicate rather heavy use of acid (hydrochloric acid) or products handled at a high temperature. Ensuring the transparency of the panels is a part of regular maintenance for the enclosure.

CLEANING THE ENCLOSURE

The enclosure must be cleaned on a regular basis.

This can be done using several different methods:

- Using soapy water then rinsing with clean water and drying with a smooth, non-abrasive paper towel.
- Using a commercial pH neutralizer and drying with a smooth, non-abrasive paper towel.
- Using a commercial window cleaner.

Metal-plated parts:

These parts should be inspected and must be free from any sign of corrosion. Verify that there is no stagnant water in the retention tray. Clean the retention tray if necessary.

WARNING

WARNING

The E.S.P.[®] program (Erlab Safety Program) was established to guarantee your safety. We remind you of the need to verify the safety parameters before using the device for the first time and in the event of a new application.

Substances that are carcinogenic, mutagenic, or toxic for reproduction (CMRs) may be handled underneath a fume hood according to the French Labor Code. However, the nature of these chemicals requires special attention from their users.

The device is not designed to be used in an explosive environment.

The filters delivered with this device must be removed from their packaging, positioned correctly, and must be suitable for the type of chemicals being handled in order to guarantee user safety.

Erlab[®] recommends testing for filter saturation on a regular basis.

Although the Timer[™] alarm can be deactivated, Erlab[®] recommends keeping the alarm activated. This alarm is triggered after every sixty hours of operation.

Erlab[®] recommends regularly monitoring the air face velocity according to the values set forth by AFNOR NF X 15-211: 2009 standard.

The quantities of the chemicals handled in the enclosure should not be greater than those listed in the guide of approved chemicals (the Chemical Listing).

AFNOR NF X 15-211: 2009 standard specifies that only those chemicals with a TLV (Threshold Limit Value) should be handled in the enclosure. These chemicals must also be included in the exhaustive list of authorized substances provided by erlab[®]. In accordance with NF X 15-211: 2009 standard, it is only possible to carry out operations that can be immediately stopped in a class 2 enclosure. The fume hood's filter must also be replaced as soon as the presence of a chemical is detected downstream from the filter.

New filters must be stored in their packaging, laid flat, in a dry place.

Erlab[®] recommends keeping a log specific to this device that includes the chemicals handled, the frequency with which these chemicals are handled, as well as any operations related to device maintenance.

**Feel free to contact our company at any time
with any questions related to the set-up,
maintenance, or use of this device.
Our team is at your disposal, providing you with clear and
detailed responses in a timely manner.**



www.erlab.com

EUROPE

FRANCE :



erlab D.F.S. S.A.S.
Parc d'Affaires des Portes BP 403
27104 Val de Reuil Cedex
Tel. : +33 (0)2 32 09 55 80
Fax. : +33 (0)2 32 09 55 90
E-Mail : Ventess@erlab.net

U.K. :



erlab D.F.S. S.A.S.
UK and Ireland representation office
Home Farm Buildings / Home Farm
Netherhampton - Salisbury - SP2 8PJ
Tel. : +44 (0)1722 341 940
Fax. : +44 (0)1722 341 950
E-Mail : SalesUK@erlab.net

GERMANY :



erlab D.F.S. S.A.S.
Vertretungsbüro Deutschland
Siegburger Straße 215
50679 Köln
Tel. : (0)800 330 47 31
Fax. : (0)800 330 47 32
E-Mail : Verkauf@erlab.net



ITALY :



erlab D.F.S. S.A.S.
Ufficio di rappresentanza in Italia
Via Leone XIII, 10 - 20145 Milano
Tel. : +39 (0)2 89 00 771
Fax. : +39 (0)2 72 097 812
E-Mail : Vendite@erlab.net

SPAIN :



erlab S.L.
Pol. Ind. Sur
Passaje Newton 3A
08754 El Papiol-Barcelona
Tel. : +34 93 673 24 74
Fax. : +34 93 673 24 76
E-Mail : Ventess@erlab.net

NORTH AMERICA

U.S.A. :



erlab inc.
388 Newburyport Turnpike
Rowley, MA 01969
Tel. : +1 (978) 948-2216
Fax. : +1 (978) 948-3354
E-mail : captainsales@erlab.com

ASIA

CHINA :



Kunshan erlab D.F.S. co Ltd.
100 Liu Shi Jing, road
Kunshan Development Zone
Jiangsu Province
Penglang - P.R. China 215333
Tel. : +86 (0) 512 5781 4085
Fax. : +86 (0) 512 5781 4082
E-mail : Sales.china@erlab.com.cn

MALAYSIA :



erlab asia sdn bhd
25 Jalan Firma - /1 Kawasan Perindustrian
Tebrau - 81 100 Johor Bahru, Johor State
Tel. : +60 (0)7 3 555 724
Fax. : +60 (0)7 3 552 810
E-Mail : Erlab@tm.net.my